PHOTOBIOMODULATION THERAPY

MiniPort Home

User Manual

Meditech International Inc
 411 Horner Ave Unit 1
 Toronto, ON M8W 4W3
 Canada

Toll Free: 1-844-770-0177 (North America Only)
 Phone: 1-416-251-1055 or 1-888-557-4004 Ext: 117 & 137
 Email: service@bioflexlaser.com
 Website: www.bioflexlaser.com

Family: BIOFLEX□ MiniPort Home
 Models: Main Controller Unit I Home (MCU I-H), DUO+ 180 (LS-B 2250+), DUO+ 240 (LS-B 3000+), DUOFLEX+ 120 (LS-B 3270), DUOFLEX+ 180 (LS-B 4890)
 Date: August 2025

Regulatory Representatives		
Region	Representative	Contact
Europe	EU REP Authorized Representative	CMC Medical Devices & Drugs SL C/ Horacio Lengo N18 29006, Málaga, Spain
United Kingdom	UK REP UK Responsible Person	CMC Medical Devices Ltd Office 32, 19-21 Crawford St London, W1H 1PJ United Kingdom
Australia	Australian Sponsor	RQSolutions Pty Ltd 44 Del Monte Pl Copacabana NSW 2251 Australia

Table of Contents

Intended Use and Indications for Use.....	4
Target Population	4
Contraindications	4
Side Effects.....	4
Important Information	4
LED Array Cooling	6
Essential Performance	6
System Components.....	6
System Setup	8
System Features.....	9
System Operation, Indicators, and Controls.....	10
Application of Treatment.....	12
BIOFLEX Practitioner+	13
Calibration, Maintenance, and Service	22
System Troubleshooting	24
System Specifications	25
Symbols Glossary	26
Safety Standards.....	27
Electromagnetic Compatibility (EMC)	28
Contact Us.....	32

Intended Use and Indications for Use

The BIOFLEX® MiniPort Home is used by trained healthcare professionals as well as general public and is indicated for use to emit energy in the red and infrared spectrum to provide topical heating for the purpose of elevating tissue temperature for temporary relief of minor muscle and joint pain, arthritis, muscle spasm, relieving stiffness, promoting relaxation of muscle tissue and temporary increase in local blood circulation.

Target Population

Adults (18+)

Contraindications

The use of BIOFLEX® systems is contraindicated for patients with lupus and other diseases which may be stimulated by light, history of keloids, recent excessive sun exposure and/or artificial tanning at the treatment area, recent patients using photosensitive medications or herbs, and patients who are pregnant. For additional advice consult your personal physician or consultants at Meditech.

Side Effects

The following may be experienced: pain; bruising and swelling; redness and inflammation; herpes simplex outbreaks and bacterial infections; temporary skin lightening or darkening; darkening/lightening of tattoos; freckle loss or lightening of moles; permanent skin pigment changes; exacerbation of preexisting skin conditions; and hair loss at the treatment area.

Important Information

- It is appropriate to read and understand the user manual. Failure to follow the instructions may result in complications and void the product warranty.
- Ensure that the system is electrically rated to operate with the alternating current power available in your location. Call BIOFLEX® Service for additional assistance.
- Connecting electrical equipment to a power bar effectively leads to creating a Medical Electrical System and the result can be a reduced level of safety (see IEC 60601-1).
- Use of accessories, detachable parts, and materials other than those provided with the system is not allowed and may be unsafe. Do not connect equipment other than what is provided with the system.
- Warning: Do not attempt to modify or open any part of the system. No modification of this equipment is allowed.

- The system is intended to be used externally over the skin only. Do not treat over compromised skin, mucous membranes, or over the eyes.
- Do not use the system over the eyes. It is inappropriate to do this as it may irritate or even damage the eyes.
- Avoid looking directly into the LED beam or staring at its bright reflections.
- Use dark glasses if your eyes are sensitive to the light emanating from the LED array or while treating upper body areas.
- Keep the system away from radiators and heat sources. Do not block any ports or cooling vents.
- Prolonged exposure to direct sunlight or excessive exposure to lint or dust is not recommended and may accelerate deterioration of the system components.
- The device must be used in a dry environment and be kept away from moisture at all times.
- The system is not waterproof. Do not allow it to become wet, do not rinse under water, or use while bathing.
- The system will produce an audible signal when powered on. Ensure the audible signal is heard before operating the system.
- Do not cover the LED array during treatment. The top surface removes excess heat into the surrounding environment.
- All system cables should be periodically inspected for signs of deterioration. Do not use the system if any cables are damaged.
- Do not exert excessive pressure or twisting motions on the LED array as this may damage the system circuitry.
- Do not pull or wrap cables around the LED array too tightly.
- Strangulation Hazard: Long cords and straps. Do not wrap cords or straps around the neck.
- Do not use the system while operating a vehicle or machinery.
- Do not use the system as a passenger in a vehicle. The bright light may distract the driver.
- Do not leave the device unattended while in operational mode.
- Keep out of reach of pets, pests, and children. Do not let young children use the system unattended.
- The system is designed for a service life of 10 years assuming that recommendations for annual calibration and maintenance are followed.
- Follow applicable local medical and electronic equipment disposal regulations and acts, including the Waste Electrical and Electronic Equipment Directive (WEEE Directive), if applicable.
- Do not run any other applications on the system except for BIOFLEX® Practitioner+.
- Always activate the Software Lock in BIOFLEX® Practitioner+ when the system is not in use.

LED Array Cooling

Surfaces of the LED arrays may get hot and may reach temperatures as high as 51°C at 40° ambient temperature. Caution must be exercised. This is particularly applicable to treatment steps running in infrared light. The higher the power and/or duty cycle, the higher the temperature. Hotter surfaces are also expected if treatments are run consecutively without sufficient time for the LED array to cool down.

If the administering person is different than the person receiving the treatment, they should be warned that the LED array may become uncomfortable and that they should notify the administering individual if so especially for young children, elderly patients and/or patients who are sensitive to heat. Common sense must prevail. If an individual leans back against the LED array placed on the spine, there may be lack of ventilation and overheating may occur. The LED array must be properly ventilated.

Cool the LED array for a minimum of 15 minutes at room temperature in case it gets too hot. Note that the internal temperature sensor will prevent the LED array from operating if the LED array gets very hot.

Essential Performance

According to IEC 60601-1 Ed. 3.2, Essential Performance is performance of a clinical function, other than that related to basic safety, where loss or degradation beyond the limits specified by the manufacturer results in an unacceptable risk. The Essential Performance of the BIOFLEX® MiniPort Home according to IEC 60601-1 Ed. 3.2, is to not radiate harmful levels of light to shine in the eyes of the patients, operators or bystanders.

System Components

Image	Name and Model
System components are intended to be used with the BIOFLEX® MiniPort Home only. Parts may appear visually different than as displayed.	
 The image shows a black, rectangular device with a color touchscreen display. The screen displays a human figure with red dots indicating treatment points, along with some text and icons. Below the screen are three circular buttons: a red one on the left, a white one in the center, and a green one on the right. The BIOFLEX logo is visible at the bottom of the device.	Main Controller Unit I Home MCU I-H



Image	Name and Model
	DUOFLEX+ 180 LS-B 4890
	DUO+ 240 LS-B 3000+
	Power Supply
	Power Cable
	Treatment Head Cable
	Straps

Image	Name and Model
	Carrying Case

System Setup

Instruction	Image
1. Connect the power supply into the socket in the back of the device.	
2. Connect the power cable to the power supply.	
3. Connect the power cable to a wall socket.	

Instruction	Image
4. Connect the treatment head cable into the port in the MCU I-H.	
5. Connect the treatment head cable to the LED array. The system is now ready for use.	

System Features

Port

The MCU I-H has one port available on the side. The port is used to connect any LED array compatible with the MCU I-H. The plug and socket are keyed, thus they must be correctly orientated before they will mate. Recognition of the LED array will come through the software.

Power Supply

A power supply is provided with your system and is connected to the MCU I-H to provide the power needed for operation. Only use the supplied power supply with your system. Do not connect any other power supply to your MCU I-H.

Power Cord

Use a Meditech-supplied power cord to connect the power supply to a power socket.

On/Off Button

The MCU I-H also has an on/off switch located on the front panel (illuminated). This button will turn BIOFLEX® Practitioner and the screen on and off once the power to the unit is available.

External USBs

The MCU I-H has two USB ports located on the back of the unit. The primary use of these ports is to connect a flash drive for backup. This USB port is also available to connect external accessories, such as a keyboard or mouse.

Software Installation

BIOFLEX® Practitioner+ is pre-installed on the MCU I-H. If updates are required, instructions will be provided by Meditech.

LED Array

The LED array is capable of emitting red or infrared light. The system will automatically progress from a red to an infrared step without pause for a single placement of the LED array. The system will pause automatically between infrared and red steps to allow for repositioning of the LED array.

System Operation, Indicators, and Controls

This section describes the system's audible and visual warnings for the LED arrays. It also explains how to stop the system in an emergency situation.

DUOFLEX+ 180

- Once enabled, two beeps will sound, one short and one long. The green indicator LED on the LED array Start/Stop button will start blinking, indicating that the LED array is ready to start.
- Once started by quickly double tapping the button on the LED array or tapping Play, two long beeps will sound. The green indicator LED will turn solid. The LED array will turn on, and light will be radiated from the LED array.
- Once paused by tapping the button on the LED array or tapping Pause, a long beep will sound. The green indicator LED will start blinking, and light from the LED array will stop radiating.
- Once stopped by tapping Stop, a long beep will be heard. The green indicator LED will start blinking. After confirmation to stop the prescription is made, another long beep will sound, and the green indicator LED on the will turn off.

DUO+ 240

- Once enabled, two beeps will sound, one short and one long. The green

indicator LED on the LED array will start blinking, indicating that the LED array is ready to start.


- Once started by quickly double tapping the button on the LED array or tapping Play, two long beeps will sound. Both the green and amber indicator LEDs will turn solid. The LED array will turn on and light will be radiated.
- Once paused by tapping the button on the LED array or tapping Pause, a long beep will sound and the amber indicator LED will turn off. The green indicator LED will start blinking and the light from the LED array will stop radiating.
- Once stopped by tapping Stop, a long beep will be heard and the amber indicator LED will turn off. The green indicator LED will start blinking. After confirmation to stop the prescription is made, another long beep will sound, and the green indicator LED on the will turn off.

Emergency Stop

In an emergency, press the Red button on the front of the MCU I-H. This immediately stops the emission of light from the LED array. Beeps will be heard and active treatment will be stopped. Treatment can only be resumed after confirmation on BIOFLEX® Practitioner+.

Placement of LED Arrays

The LED arrays may be secured directly to the body with the provided strap. This is done by feeding the strap through the end loop of the LED array and then mating the strap.

- Do not bend the LED array beyond its normal curvature, as this may break the LED array.
- The LED array must be secured while allowing air circulation around the top surfaces of the LED array.

Application of Treatment

1. Local: Apply the LED array to the painful area.
2. Treating a joint: When treating a joint, place the LED array over the most symptomatic area and if possible, surround the entire joint. For example, if treating hand arthritis, place the LED array on the palm and back of the hand.
3. Good contact between the LED array and skin should always be maintained. This may require gentle application of pressure keeping in mind that ventilation for proper cooling must continue. When strapping an LED array, good contact without excessive pressure is important.
4. Maintaining cleanliness: Ensure the skin is free of oil, lotions, creams, etc. before applying the LED array. LED arrays should be carefully cleaned between use.
5. Frequency of use: Initially treat once per day, but not exceeding twice a day. As symptoms improve, treatment can be reduced to every other day until symptoms resolve. For progressive conditions like osteoarthritis,

ongoing treatment 2 - 3 times a week may be helpful in managing ongoing symptoms of pain and stiffness.

Treatment Protocols and Levels


The preprogrammed protocols provide optimal parameters depending on the condition being treated. Levels of protocols reference increasing dosage of light and stimulation. It is recommended to start with Level 1 and only increase to the next level if symptoms are not improving or are no longer improving.

LED Arrays

Your BIOFLEX® MiniPort Home is provided with a DUOFLEX+ 180 by default; however, an optional DUO+ 240 is available. Both of these LED arrays provide the same dosage. Thus, any LED array can be used for any specific protocol. LED array selection depends on user preference.

Placement


The system is primarily utilized to treat localized pathologies characterized by pain and other symptoms. It may also be placed over the area of the spine innervating the area that is symptomatic. The following illustrations demonstrate typical ways to place the array over a treatment area:

BIOFLEX Practitioner+

BIOFLEX® MiniPort Home Practitioner+ (i.e. BIOFLEX® Practitioner+) is a software application developed by Meditech and designed to work with the BIOFLEX®

MiniPort Home. It is an intuitive and easy to use application. The following sections describe the use of the application and its features.

From the Home screen, choose to run a new treatment, a previous treatment, or a favourite treatment. Selecting to run a new treatment will lead to the protocol selection window. Selecting to run a previous treatment will lead to a list of previously run treatments. While running any protocol, an option exists to select that treatment as a favourite treatment. Selecting to run a favourite treatment will lead to a list of favourite treatments.

On the side bar choose to go to Protocols, Favourite, Help, Lock System, Settings, or the On/Off button (to turn the application on/off).

On the top right is a slider for adjusting the speaker volume, a button for viewing the port status and a button for additional options.

Side Bar

The side bar is located on the left of the screen. It enables the following actions:

- The Home button  can be clicked at any time to go back to the Home screen.
- The Protocols button  can be clicked at anytime to go to the Select Protocol screen. From the Protocol screen, select a protocol.

- The Favourite button  can be clicked at anytime to go to the Favourite Protocols screen. From the list of favourite protocols, choose a protocol to run.
- The Help button  can be clicked at any time to access the Help screens. The Help screens contain the information of this chapter as well as a link to the treatment tutorials, links to the frequently asked questions, and information pertaining to the software name, software version, and software region.
- The Lockup button  provides a way to lock the system in case the operator needs to step away. This eliminates the risk of an unauthorized person using the device. Once the system is locked, a password is required to unlock the system. The Lockup factory setting password is locklaser. This password can be changed by choosing the Settings button  from the side bar in the Home screen and following the steps for changing the password.
- The Settings button  can be used to change settings, such as passwords, volume, brightness, etc.
- The On/Off button  allows the user to shut down the application. Note that the device is still powered on. To fully turn the MCU I-H off, unplug the power supply located at the back of the unit.

Selecting a Protocol

A protocol can be selected after choosing to run a new treatment from the home screen or by selecting the protocols button from the side bar. The Select Protocol screen shows up.


In the Select Protocol screen:

- Choose one of the body locations (dots on the body figures) or a protocol under other conditions.
- If there is more than one protocol associated with a body location, select the most applicable protocol under Available Protocols.


- After selecting a protocol, it can be added to the Favourite Treatments list by clicking the Add Favourite button. If a protocol is already a favourite treatment the button will change to Remove Favourite.
- Make sure to select the desired level before selecting Next. There are no set rules or set levels for protocol treatments. It is recommended to start at Level 1 for treatments and evaluate their success. In this way, the levels can be raised or lowered based on treatment performance.

- Click Run Treatment to start the treatment.

Running a Treatment

After selecting a protocol, the LED array is activated (a message will come up if no LED array is connected) and the following screen is shown. The Running Treatment screen shows the progress of the treatment.

- To start the treatment click the Play button (▶).
- Treatment can be paused and resumed by clicking the Pause button (⏸) and the Play button (▶).
- Information about the energy density applied is shown as the treatment progresses.
- The pictures at the bottom of the screen are typical placements of the treatment being run; however, the exact location is determined by the location of symptoms.
- Placements can be skipped by clicking on any step in the cascading list of steps (pause a running step first). Not all placements are necessary. Actual placements are determined by the operator and based on the location and extent of the pathology being treated. Skip step(s) as necessary.

Dosage

Dosage (J/cm^2) is calculated and displayed as the treatment progresses. The user can see the energy density that is delivered during the treatment and at the end of the treatment the total dosage that was delivered.

Running a Previous or Favourite Treatment

Running a previous treatment or a favourite treatment is an option from the Home screen.

- Clicking on the Run the Previous Treatment button on the Home screen will take you to the Past Treatments screen. Select one of the treatments from the list to run.
- Similarly, clicking on the Run a Favourite Treatment button on the Home


screen will take you to the Favourite Treatments screen, which is a list of all favourite treatments. Select one of the treatments from the list to run.

- Instructions and steps after selecting a previous or favourite treatment are the same as running a selected protocol.

Port Status

A port's status can be seen at any point by clicking the Port Status button. The Port Controller Status window will show.

The Port Controller Status window shows the port status, serial number, and firmware version. To see the port data, click the Refresh Port Data button. To

reset the port, click the Reset Port button.

The Port Controller Status window will also show the model of the LED array that is connected to the port together with its output specifications. The internal temperature inside the LED array and if a treatment is running or paused, the remaining time for that treatment, are also displayed.


The port status indicates the state of the port. Some states of the treatment port can be:

- N/A: An LED array is not connected to the port.
- Idle: An LED array is connected to the port, but no protocol is selected.
- Linked: An LED array is connected to the port and a treatment protocol is selected.
- Enabled: An LED array is connected to the port and is enabled (LED array is fully powered and ready to be started).
- Running: The port is running a treatment protocol.
- Paused: The port is paused in the middle of treatment.
- Error (e.g. overheat)

If the port is in an error state, the user needs to click the Get Error Code button and proceed as instructed. The error may be cleared by clicking the Clear Error button in the Port Controller Status window.


Treatment Tutorials

Tutorial support for treatments is available from the Help button in the side bar or the Additional Options button in the Home screen. When the Tutorials screen opens, an introductory video starts (if applicable).

Similar to selecting a protocol for treatment, the Tutorial screen has body locations and other protocols to choose from. Once a protocol is selected, the corresponding tutorial will be displayed.

If there are multiple protocols associated with an anatomy (body location), the first tutorial will be shown automatically. Click the other available protocols to view their tutorial.


Settings and Configuration

Settings and configurations are available by clicking the Settings button  in the side bar. Click on one of the two tabs, User Privilege Passcode or User Preference,

to access the setting that is required.

User Privilege Passcode


The User Privilege Passcode tab is used to set or change the lock system passcode. The Set Lockup Passcode button is used to set up a password for locking the system. The default lockup factory setting passcode is locklaser.

User Preference


The User Preference tab is used to change the following settings:

- **Brightness:** Change the brightness of the screen.
- **Speaker Volume:** Change the volume of the internal speakers (used when playing tutorial videos). Note this is not the volume of the audible beeps when treatments are started, paused, or stopped. The audible beeps are a safety mechanism that cannot be changed.

Additional Options

The Additional Options button is on the top right of the screen, and it can be used to access the settings and configurations or the treatment tutorials:

Calibration, Maintenance, and Service

Warning: Do not modify this equipment without Meditech's authorization. If the BIOFLEX® MiniPort Home gives repeated errors it might need calibration. Please contact BIOFLEX® Service. To ensure that the system runs at optimum performance, it is strongly recommended that the LED arrays are checked annually by BIOFLEX® Service.

Attention: Only Meditech authorized service personnel may calibrate or service the system. Unauthorized changes to the system may be subject to the immediate termination of warranty (at Meditech's discretion). Unauthorized changes to the system-specific parameters will also cause the immediate termination of warranty (at Meditech's discretion).

To reduce the risk of electric shock, do not remove the cover from any part of the system. There are no serviceable parts inside. All system cables should be periodically inspected for signs of deterioration. If damaged, contact BIOFLEX® Service immediately to fix or replace the cable(s). On request, Meditech will make available circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist Meditech authorized service personnel to repair the system.

Cleaning

- The MCU I-H may be cleaned using a damp soft cloth or a soft disinfectant cloth.
- The LED arrays may be cleaned with alcohol wipes, a soft disinfectant cloth, a soft damp cloth, and/or alcohol and cotton swab combination. Meditech recommends using PDI Sani-Cloth Plus disinfecting wipes.
- A soft bristle brush can be used to remove particles from the LED arrays.
- It is recommended that the LED arrays be cleaned after each treatment.

Follow the instructions below to clean your system:

- To clean the MCU I-H, gently wipe surfaces as necessary.
- To clean the LED arrays, use clean wipes, cloth, and/or cotton swab. If present, remove heavy soil loads with a soft bristle brush prior to cleaning. Unfold a clean wipe (or cloth) and thoroughly wipe the surfaces of the LED arrays. Make sure to thoroughly wipe the LED surfaces. Make sure to wet all surfaces and not backtrack after wiping. If necessary, use a cotton swab and alcohol to get into smaller surfaces. Allow LED arrays to air dry prior to use.

Note: Do not spray or immerse the LED arrays with water/disinfectants/cleaning solutions or wrap them with disinfectant-soaked towels. Excessive disinfecting or cleaning solutions can penetrate the equipment housings and will likely cause corrosion of the electronic circuitry, causing LED array failure.

System Damage

If part or all of the system is dropped or damaged, check it to make sure it still functions properly. To troubleshoot a damaged system, perform the following steps:

1. Check all power and interface cable connections. Make sure all cables are secure and properly connected.
2. Turn on all components of the system. If the system does not work, contact BIOFLEX® Service.

Liquids from spills, splashes, and excessive humidity can also damage your system. If the system becomes wet, perform the following steps:

1. Immediately unplug all system components from the electrical outlet.
2. Contact BIOFLEX® Service.

System Troubleshooting

Contact BIOFLEX® Service if you are unable to solve the problem with the information below:

Error Code	Problem	Solution
E10, E11, E12, E13, E14	This port controller board has failed.	1. Disconnect and reconnect the treatment head cable. 2. Replace the treatment head cable. 3. Contact BIOFLEX® Service.
S01, S02	LED array cannot switch light source or color.	Contact BIOFLEX® Service.
S04	Unable to turn on the treatment head on the port.	1. Disconnect and reconnect the treatment head cable. 2. Replace the treatment head cable. 3. Contact BIOFLEX® Service.
S06	Temperature sensor for the LED array connected to the port has failed.	1. Wait until the LED array cools to a safe temperature. 2. Resume treatment.
S07	LED array connected to the port has a factory calibration error.	Contact BIOFLEX® Service.

Error Code	Problem	Solution
S12	The LED array ran into a communication error.	1. Reset the system by unplugging the power supply, waiting 10 sec, and then replugging it in. 2. Disconnect and reconnect the treatment head cable. 3. Connect the LED array to a known working system. 4. Replace the treatment head cable. 5. Contact BIOFLEX® Service.

System Specifications

Electrical Specifications	
Rated Voltage	100 to 240 V
Rated Frequency	50 to 60 Hz
Operating Conditions (Note: Also applies to LED arrays.)	
Temperature	+5 to +40°C
Relative Humidity	15 to 93% non-condensing
Atmospheric Pressure	700 to 1,060 hPa
Storage Conditions (Note: Also applies to LED arrays.)	
Temperature	-25 to +70°C
Relative Humidity	Up to 93% non-condensing
Functional Specifications	
Waveform	Square
Percent of Total Power	0 to 100% (1% steps)
Duration	0 to 3,600 sec (60 min)
Frequency (Modulation Mode)	1 to 10,000 Hz (5 µs)
Duty Cycle	1 to 99% (1% steps, 5 µs)
Classification	

The system is classified as a medical device. In accordance to IEC 60601-1:2005+A2:2020 (3.2 Edition):
 Class II equipment, according to the type of protection against electric shock;
 Type BF applied parts, according to the degree of protection against electric shock;
 Classified IP21, according to the degree of protection against harmful ingress of water or particulate matter;
 Continuous operation, according to the mode of operation;
 None of the parts or the system is sterilized.

LED Array Specifications

Parameter	DUOFLEX+ 180	DUO+ 240
Model	LS-B 4890	LS-B 3000+
LED Maximum Optical Power (mW)	R: 9.3 IR: 18.5	R: 4.2 IR: 8.5
Array Optical Power (mW)	R: 1630 IR: 3260	R: 1000 IR: 2000
Effective Surface Area (cm ²)	163	100
Maximum Power Density (mW/cm ²)	R: 10.0 IR: 20.0	
Number of LEDs	176	240
Medium	R: AlGaInP IR: AlGaAs	
Wavelength (nm)	R: 660 IR: 840	
Spectral Width (nm)	R: 32 IR: 80	

Symbols Glossary

Symbol	Title and Description	Standard
	Refer to instruction manual/booklet	ISO 7010
	Manufacturer Indicates the medical device manufacturer.	ISO 15223-1
	Date of manufacture Indicates the date when the medical device was manufactured.	ISO 15223-1
	Authorized representative in the European Community Indicates the Authorized representative in the European Community.	ISO 15223-1
	Authorized representative in the UK Indicates the Authorized representative in the United Kingdom.	-
	Medical device Indicates the item is a medical device	ISO 15223-1

Symbol	Title and Description	Standard
	Serial number Indicates the manufacturer's serial number so that a specific medical device can be identified.	ISO 15223-1
	Warning; Hot surface	ISO 7010
	Waste Electrical and Electronic Equipment This electrical and electronic equipment is subject to separate waste collection.	EN 50419
	Type BF applied part To identify a type BF applied part complying with IEC 60601-1.	IEC 60417
	ETL listed mark Indicates that the product has been tested by Intertek and found in compliance with accepted national standards.	-
IP21	Degrees of protection provided by enclosures (IP Code)	IEC 60529

Safety Standards

This product is medical equipment and is in compliance with the following standards:

General Basic Safety:

- IEC 60601-1:2005+A2:2020 (3.2 Edition). Canada and USA national deviations:
 - CAN/CSA-C22.2 No. 60601-1:2014+A2:2022 (3.2 Edition)
 - ANSI/AAMI ES60601-1:2005+A2:2021 (3.2 Edition)
- IEC 60601-1-11:2015+A1:2020 (2.1 Edition)
- IEC 60601-1-6:2010+A1:2013+A2:2020 (3.2 Edition)
- 21 CFR 1040 - Performance Standards for Light-Emitting Products

EMC:

- IEC 60601-1-2:2014+A1:2020 (4.1 Edition)

LED Safety:

- IEC 60601-2-57:2011 (1st Edition)
- IEC 62471:2006 (1st Edition) including European deviations

Electromagnetic Compatibility (EMC)

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual. Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

List of Cables, Transducers, and Accessories			
Cable	Length	Connected From	Connected To
Power Cord, CT402.019	1.83 m	Power Supply	Power Outlet

Warning: The use of accessories and cables other than those specified in this manual, with the exception of accessories and cables qualified and sold by Meditech International Inc may result in increased emissions or decreased immunity of the equipment and may cause the system to be non-compliant with the requirements of IEC 60601-1-2:2014+A1:2020.

Warning: The BIOFLEX® MiniPort Home should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment should be observed to verify normal operation in the configuration in which it will be used.

Example Emission Class and Group Compliance		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
The BIOFLEX® MiniPort Home is intended for use in the electromagnetic environment specified below. The customer or user of the system should assure that it is used in such an environment.		
RF Emissions CISPR 11	Group 1	This system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Example Emission Class and Group Compliance

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Class B	The emissions characteristics of this equipment make it suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

Immunity and Test Level Compliance

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
The BIOFLEX® MiniPort Home is intended for use in the electromagnetic environment specified below. The customer or user of the system should assure that it is used in such an environment.			
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV Contact ± 15 kV Air	± 8 kV Contact ± 15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/ Burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line to line	± 0.5 kV, ± 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.

Immunity and Test Level Compliance			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Voltage Dips, Short Interruptions, and Voltage Variations on Power Supply Input Lines IEC 61000-4-11	0% U_T (100% dip in U_T) for 0.5 cycle 0% U_T (100% dip in U_T) for 1 cycles 70% U_T (30% dip in U_T) for 25/30 cycles 0% U_T (100% dip in U_T) for 5 sec	0% U_T (100% dip in U_T) for 0.5 cycle 0% U_T (100% dip in U_T) for 1 cycles 70% U_T (30% dip in U_T) for 25/30 cycles 0% U_T (100% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of this system requires continued operation during power mains interruptions, it is recommended that the system be powered from an uninterruptible power supply or a battery.
Power Frequency Magnetic Field (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM/Amateur Radio bands inside 150 kHz to 80 MHz 3 V/m 80 MHz to 2.7 GHz RF communication equipment inside 80 MHz to 6 GHz	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM/Amateur Radio bands inside 150 kHz to 80 MHz 3 V/m 80 MHz to 2.7 GHz RF communication equipment inside 80 MHz to 6 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the system including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .

Immunity and Test Level Compliance			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Note 1: U_1 is the AC mains voltage prior to application of the test level. Note 2: At 80 MHz and 800 MHz, the higher frequency range applies. Note 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the system. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the BIOFLEX® MiniPort Home			
The BIOFLEX® MiniPort Home is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.24
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of these systems, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Contact Us

Please contact BIOFLEX® Service and the competent authority of the EU member state in which you or the patient is established if any serious incident has occurred in relation to the device.

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