



MiniPort Professional

The latest edition to the BIOFLEX family. This single port system offers a lightweight and portable laser therapy experience. Perfect for the clinician on the go.

Includes:

- + Touch screen monitor
- + BIOFLEX DUO+240 Treatment Array
- + BIOFLEX Infrared Laser Probe
- + One set of safety glasses

Features:

- + State of the art patient database allows monitoring and tracking patient progress on four screens simultaneously
- + Intuitive touchscreen that continuously monitors treatments throughout the entire session
- + Manual stop button for increased safety
- + Surprisingly compact for a full feature system

Technical Specifications

Function	Specification
Mode of operation	CW, Pulsed (Modulated)
Waveform	Square
Percent of total power	0-100% (1% steps)
Duration	0-3,600 seconds (60:00 minutes)
Frequency	1-10,000 Hz (5 μ s resolution)
Duty Cycle	10-90% (10% steps) (5 μ s resolution)

Patents + Approvals

Compliance with the requirements to bear TÜV SÜD marks.

BIOFLEX Systems are manufactured under the ISO 13485:2016 quality system.

BIOFLEX MiniPort Systems are authorized for sale by Health Canada and registered with the FDA as a medical device.

BIOFLEX systems are approved for marketing in Canada under a Health Canada Medical Device License.

Meditech holds patents worldwide, including but not limited to North America, Europe and Australia