



MultiPort Professional

The BIOFLEX® MultiPort System can operate up to four treatment arrays or laser probes simultaneously. Continuous monitoring is displayed on all four screens throughout the entire treatment.

Includes:

- + Touch screen monitor
- + Two BIOFLEX DUO+240 Treatment Arrays
- + Two laser probes of your choice
- + Two sets of safety glasses

Features:

- + Treat up to 4 patients simultaneously using any combination of flexible arrays and laser probes
- + Intuitive touchscreen that continuously monitors treatments throughout the entire session
- + Manual stop button for increased safety
- + Our most advanced, full featured Laser Therapy System

Technical Specifications

Function	Specification
Mode of operation	CW, Pulsed (Modulated)
Waveform	Square
Percent of total power	0-100% (1% steps)
Duration	0-3,600 seconds (60:00 minutes)
Frequency	1-10,000 Hz (5 μ s resolution)
Duty Cycle	10-90% (10% steps) (5 μ s resolution)

Patents + Approvals

Compliance with the requirements to bear TÜV SÜD marks.

BIOFLEX Systems are manufactured under the ISO 13485:2016 quality system.

BIOFLEX MultiPort Systems are CE marked and registered with the FDA in the United States.

BIOFLEX systems are approved for marketing in Canada under a Health Canada Medical Device License.

Meditech holds patents worldwide, including but not limited to North America, Europe and Australia.